

Anti-HOXA1 antibody (Internal) (STJ72152)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/WB-Trf
Host / Source	Goat
Reactivity	Human/Mouse/Rat/Dog/Cow

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:8000. WB: In transfected HEK293 transiently expressing Human HOXA1 (MYC and DYKDDDDK tagged) , a band of approx. 19kDa is observed. No bands are observed in mock-transfected HEK293 and the same b
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3198
Gene Symbol	HOXA1
UniProt ID	HXA1_HUMAN
Immunogen Region	Internal
Immunogen Sequence	QPATYQTSGN
Specificity	This antibody is expected to recognize reported isoforms a and b (NP_005513.1; NP_705873.2).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



HEK293 lysate (10ug protein in RIPA buffer) overexpressing Human HOXA1 (RC222721) with C-terminal MYC tag probed with STJ72152 (1ug/ml) in Lane A and probed with anti-MYC Tag (1/1000) in lane C. Mock-transfected HEK293 probed with STJ72152 (1mg/ml) in Lane B. Primary incubations were for 1 hour. Detected by chemiluminescence.