

Anti-SLC4A8 antibody (Internal) (STJ72197)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	Pep-ELISA
Host / Source	Goat
Reactivity	Human/Mouse/Rat/Dog/Cow
PRODUCT PROPERTIES	
Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:1000. WB: Preliminary experiments in Human Brain (Cerebellum) , Mouse Brain and Rat Brain lysates, and in lysates of Human Breast Cancer, Human Kidney and Human Lymph Node gave no specific signal
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	9498
Gene Symbol	SLC4A8
UniProt ID	S4A8_HUMAN
Immunogen Region	Internal
Immunogen Sequence	TNLEVKNGVNCEH
Specificity	This antibody is expected to recognize isoform a and b (NP_001035049.1; NP_004849.2).
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.