

Anti-FANCF antibody (217-228) (STJ72279)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/WB-Trf
Host / Source	Goat
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:16000. WB: In transfected HEK293 transiently expressing full-length Human FANCF (myc and DYKDDDDK tagged) , a band of approx. 40kDa was observed. No bands were observed in mock-transfected HEK293
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2188
Gene Symbol	FANCF
UniProt ID	FANCF_HUMAN
Immunogen Region	217-228
Immunogen Sequence	RPQEELEPGIHK

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



HEK293 lysate (10ug protein in RIPA buffer) overexpressing Human FANCF with C-terminal MYC tag probed with STJ72279 (1ug/ml) in Lane A and probed with NA anti-MYC Tag (1/1000) in lane C. Mock-transfected HEK293 probed with STJ72279 (1mg/ml) in Lane B. Primary incubations were for 1 hour. Detected by NA chemiluminescence.