

Anti-BIRC6/Apollon antibody (4724-4735) (STJ72374)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA
Host / Source	Goat
Reactivity	Human/Mouse/Rat/Cow

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:64000. WB: Preliminary experiments gave bands at approx 90kDa and 70kDa in lysates of cell lines A431 and HeLa after 0.1µg/ml antibody staining. Please note that currently we cannot find an expla
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	57448
Gene Symbol	BIRC6
UniProt ID	BIRC6_HUMAN
Immunogen Region	4724-4735
Immunogen Sequence	SREYDGNIRQAT

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---