

Anti-NRAP antibody (Internal) (STJ72382)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA
Host / Source	Goat
Reactivity	Human/Mouse/Rat/Dog/Cow/Pig

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:4000. WB: Preliminary experiments in Human Heart, Skeletal Muscle and Uterus lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any f
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4892
Gene Symbol	NRAP
UniProt ID	NRAP_HUMAN
Immunogen Region	Internal
Immunogen Sequence	EDTKSQCHVSLD
Specificity	This antibody is expected to recognize all reported isoforms (NP_006166.3; NP_932326.2).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---