

Anti-BLK antibody (N-Term) (STJ72502)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/WB/IHC
Host / Source	Goat
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. WB: Approx 50kDa band observed in Human Thymus, Spleen and Lymph Node lysates (calculated MW of 57.7kDa according to NP_001706.2). Recommended concentration: 0.1-0.3µg/ml. Primary incubat
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	640
Gene Symbol	BLK
UniProt ID	BLK_HUMAN
Immunogen Region	N-Term
Immunogen Sequence	QWSPLKVSAQDKD

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



STJ72502 (0.1µg/ml) staining of Human Thymus lysate (35µg protein in RIPA buffer). Detected by chemiluminescence.