

**Anti-DEFB107A/DEFB107B antibody (Internal) (STJ72659)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	Pep-ELISA/IHC
Host / Source	Goat
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. WB: Not yet tested-our routinely used western blotting protocol does not allow detection of proteins as small as the calculated size of 7.85kDa according to NP_001032757.2. IHC: Paraff
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

**TARGET INFORMATION**

Immunogen Region	Internal
Immunogen Sequence	HRALISKRMEGHCE

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---