

Anti-Acacb-Mouse antibody (Internal) (STJ72669)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	Pep-ELISA/IHC
	Host / Source	Goat
	Reactivity	Mouse/Rat/Pig
PRODUCT PROPERTIES		
	Clonality	Polyclonal
	Concentration	0.5 mg/mL
	Conjugation	Unconjugated
	Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
	Dilution Range	Peptide ELISA: antibody detection limit dilution 1:16000. WB: Not yet tested-our routinely used western blotting protocol does not allow detection of proteins as large as the calculated size of 277kDa according to NP_001084.3. Therefore we cannot
	Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
	Isotype	IgG
	Storage Instruction	Store at-20°C on receipt and minimise freeze-thaw cycles.
TARGET INFORMATION		
	Immunogen Region	Internal
	Immunogen Sequence	TLNVQESGRVKYIKR
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.