

Anti-DBH antibody (437-448) (STJ72732)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/WB/IHC
Host / Source	Goat
Reactivity	Human/Mouse/Rat/Pig/Cow

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. WB: Approx 70kDa band observed in lysates of cell line Jurkat (calculated MW of 69.1kDa according to NP_000778.3). Recommended concentration: 1-3µg/ml. Primary incubation was 1 hour. Prel
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1621
Gene Symbol	DBH
UniProt ID	DOPO_HUMAN
Immunogen Region	437-448
Immunogen Sequence	DNHYSPHFQEIR

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



STJ72732 (2µg/ml) staining of Jurkat lysate (35µg protein in RIPA buffer). Detected by chemiluminescence.