

Anti-Praja-1/RNF70 antibody (285-298) (STJ72775)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA
Host / Source	Goat
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:4000. WB: Preliminary experiments gave an approx 30kDa band in Human Brain and Mouse Brain lysates after 1µg/ml antibody staining. Please note that currently we cannot find an explanation in the
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	64219
Gene Symbol	PJA1
UniProt ID	PJA1_HUMAN
Immunogen Region	285-298
Immunogen Sequence	HSTSRWRDTANDNE
Specificity	This antibody is expected to recognize all reported isoforms (NP_660095.1; NP_001027568.1; NP_071763.2).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---