

Anti-WDR9/BRWD1 antibody (2137-2151) (STJ72789)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA
Host / Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. WB: Preliminary experiments gave an approx 100kDa band in Human Heart and Skeletal Muscle lysates after 0.3µg/ml antibody staining. Please note that currently we cannot find an explanation
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54014
Gene Symbol	BRWD1
UniProt ID	BRWD1_HUMAN
Immunogen Region	2137-2151
Immunogen Sequence	NVKISETTGNSKFRP
Specificity	This antibody is expected to recognize reported isoforms A and B (NP_061836.2; NP_387505.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---