

Anti-E2F1 antibody (480-493) (STJ72874)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/IF
Host / Source	Goat
Reactivity	D.melanogaster

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. IF: This antibody shows identical pattern in the Drosophila eye-antennal imaginal disc as was observed before with a different antibody (Moon et al, Mol Cell Biol. 2006 Oct;26 (20) :7601-15
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	480-493
Immunogen Sequence	TKGLFQTPYRSARN
Specificity	Reported isoforms represent identical protein: NP_524437.2; NP_732646.1; NP_732647.1

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



STJ72874 (0.5 µg/ml) staining of the Drosophila Eye-Antennal Imaginal Disc. Data obtained from Dr. Nam-Sung Moon, McGill University, Montreal, Canada.