

Anti-THNSL2 antibody (Internal) (STJ72925)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/IHC
Host / Source	Goat
Reactivity	Human/Dog/Cow

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. WB: Preliminary experiments gave an approx 35kDa band in Human Colorectal Cancer lysates after 1µg/ml antibody staining. Please note that currently we cannot find an explanation in the li
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55258
Gene Symbol	THNSL2
UniProt ID	THNS2_HUMAN
Immunogen Region	Internal
Immunogen Sequence	SRFRHREVVHLSR
Specificity	This antibody is expected to recognize isoform 1 (NP_060741.3) and isoform 2 (NP_001231605.1) only.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---