

Anti-TNNI2 antibody (N-Term) (STJ72983)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/WB
Host / Source	Goat
Reactivity	Human/Mouse/Rat/Dog

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. WB: Approx. 25kDa band observed in Human, Mouse and Rat Skeletal Muscle lysates (calculated MW of 21.3kDa according to NP_003273.1). Recommended concentration: 0.03-0.1µg/ml. Primary incu
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7136
Gene Symbol	TNNI2
UniProt ID	TNNI2_HUMAN
Immunogen Region	N-Term
Immunogen Sequence	KRNRAITARRQHLKS
Specificity	This antibody is expected to recognize both reported isoforms (NP_003273.1; NP_001139313.1). Reported variants represent identical protein: NP_001139301.1, NP_003273.1.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



STJ72983 (0.03µg/ml) staining of Human (A) , Mouse (B) and Rat (C) Skeletal Muscle lysate (35µg protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.