

Anti-AHSG antibody (Internal) (STJ73038)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/WB
Host / Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. WB: Approx 60kDa band observed in Human albumin-depleted Serum lysates (calculated MW of 39.3kDa according to NP_001613.2). This molecular weight is routinely observed by other sources. R
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	197
Gene Symbol	AHSG
UniProt ID	FETUA_HUMAN
Immunogen Region	Internal
Immunogen Sequence	QPEGANEAVPTP

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



STJ73038 (0.5 µg/ml) staining of Human albumin-depleted Serum (35 µg protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.