

Anti-CD96 antibody (C-Term) (STJ73086)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	Pep-ELISA
Host / Source	Goat
Reactivity	Human
PRODUCT PROPERTIES	
Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. WB: Preliminary experiments gave an approx 100+75kDa bands in Jurkat lysates after 0.3µg/ml antibody staining. Please note that currently we cannot find an explanation in the literature f
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	10225
Gene Symbol	CD96
UniProt ID	TACT_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	KYTCIQEPNESD
Specificity	This antibody is expected to recognize both reported isoforms (NP_937839.1; NP_005807.1).
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.