

Anti-Nrac-Mouse antibody (17-31) (STJ73117)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	Pep-ELISA
	Host / Source	Goat
	Reactivity	Mouse
PRODUCT PROPERTIES		
	Clonality	Polyclonal
	Concentration	0.5 mg/mL
	Conjugation	Unconjugated
	Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
	Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. WB: Preliminary experiments in Mouse Adipose and Heart lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from p
	Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
	Isotype	IgG
	Storage Instruction	Store at-20°C on receipt and minimise freeze-thaw cycles.
TARGET INFORMATION		
	Immunogen Region	17-31
	Immunogen Sequence	CHPTRENEGTTGRSQ
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.