

Anti-SRRT/ARS2 antibody (C-Term) (STJ73145)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/WB/IHC
Host / Source	Goat
Reactivity	Human/Mouse/Rat/Dog/Pig/Cow

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. WB: Approx 140kDa band observed in lysates of cell line Daudi and of Mouse Spleen, and approx 150kDa band observed in lysates of Rat Testis (calculated MW of 101kDa according to Human NP_
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION


Gene ID	51593
Gene Symbol	SRRT
UniProt ID	SRRT_HUMAN
Immunogen Region	C-Term
Immunogen Sequence	DPRAIVEYRDLAPD
Specificity	This antibody is expected to recognize all reported isoform (NP_056992.4; NP_001122324.1; NP_001122325.1; NP_001122326.1).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



STJ73145 (0.1 µg/ml) staining of Daudi (A) and Mouse Spleen (B) lysates (35 µg protein in RIPA buffer). Detected by chemiluminescence.



STJ73145 (0.1 µg/ml) staining of Rat Testis lysate (35 µg protein in RIPA buffer). Detected by chemiluminescence.