

Anti-U-PAR/CD87 antibody (Internal) (STJ73158)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	Pep-ELISA/WB
Host / Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. WB: Approx 60kDa band observed in lysates of cell line U251 (calculated MW of 37.0kDa according to NP_002650.1). The observed molecular weight corresponds to the glycosylated form. Recomm
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5329
Gene Symbol	PLAUR
UniProt ID	UPAR_HUMAN
Immunogen Region	Internal
Immunogen Sequence	RTGLKITSLTE
Specificity	This antibody is expected to recognize all reported isoforms (NP_002650.1; NP_001005376.1; NP_001005377.1).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



STJ73158 (1 µg/ml) staining of U251 lysate (35 µg protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.