

## Anti-ANXA2 antibody (88-137 aa) (STJ91601)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/IF/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 IHC 1:100-1:300 ELISA 1:10000 IF 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Molecular Weight	Observed: 38kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION


Gene ID	<a href="#">302</a>
Gene Symbol	<a href="#">ANXA2</a>
UniProt ID	<a href="#">ANXA2_HUMAN</a>
Immunogen Region	88-137 aa
Specificity	Annexin II Polyclonal Antibody detects endogenous levels of Annexin II protein.

### ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Western blot analysis of lysate from HeLa cells, using Annexin II antibody.



Western blot analysis of HeLa cells using Annexin II Polyclonal Antibody diluted at 1: 2000



Immunohistochemistry analysis of Annexin II antibody in paraffin-embedded human brain tissue.



Western blot analysis of various cells using Annexin II Polyclonal Antibody diluted at 1: 2000