

Anti-BMP8B antibody (261-310 aa) (STJ91871)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 ELISA 1:40000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Molecular Weight	Observed: 45kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	656
Gene Symbol	BMP8B
UniProt ID	BMP8B_HUMAN
Immunogen Region	261-310 aa
Specificity	BMP-8B Polyclonal Antibody detects endogenous levels of BMP-8B protein.

ADDITIONAL INFORMATION


Note: STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of various cells using BMP-8B Polyclonal Antibody diluted at 1:2000. The lane on the right was blocked with STJ91871.



Western blot analysis of lysates from 293 and A549 cells, using BMP8B Antibody. The lane on the right was blocked with STJ91871.



Immunofluorescence analysis of SiHa cells. 1, Primary antibody was diluted at 1:100 (4°C overnight). 2, Goat Anti Rabbit STJ91871.