

Anti-BNIP2 antibody (73-122 aa) (STJ91876)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/IF/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 IHC 1:100-1:300 ELISA 1:40000 IF 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Molecular Weight	Observed: 36kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	663
Gene Symbol	BNIP2
UniProt ID	BNIP2_HUMAN
Immunogen Region	73-122 aa
Specificity	BNIP-2 Polyclonal Antibody detects endogenous levels of BNIP-2 protein.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Western blot analysis of lysate from Jurkat cells, using BNIP-2 antibody.



Western blot analysis of K562 cells using BNIP-2 Polyclonal Antibody



Immunohistochemistry analysis of BNIP-2 antibody in paraffin-embedded human breast carcinoma tissue.



Western blot analysis of various cells using BNIP-2 Polyclonal Antibody