

Anti-HMG20B antibody (1-50 aa) (STJ91880)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/IF/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 IHC 1:100-1:300 IF 1:200-1:1000 ELISA 1:20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Molecular Weight	Observed: 36kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	10362
Gene Symbol	HMG20B
UniProt ID	HM20B_HUMAN
Immunogen Region	1-50 aa
Specificity	BRAF35 Polyclonal Antibody detects endogenous levels of BRAF35 protein.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Western blot analysis of various cells using BRAF35 Polyclonal Antibody diluted at 1:2000, cells nucleus extracted by STJ91880



Western blot analysis of AD293 cells using BRAF35 polyclonal antibody diluted at 1:2000, cells nucleus extracted by STJ91880



Western blot analysis of lysates from Jurkat cells, using HMG20B Antibody. The lane on the right was blocked with the STJ91880



Immunohistochemical analysis of paraffin-embedded human spleen. 1, Antibody was diluted at 1:200 (4°C overnight). 2, STJ91880