

Anti-CLCNKA antibody (581-630 aa) (STJ92328)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 ELISA 1:40000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Molecular Weight	Observed: 75kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	1187
Gene Symbol	CLCNKA
UniProt ID	CLCKA_HUMAN
Immunogen Region	581-630 aa
Specificity	CLC-KA Polyclonal Antibody detects endogenous levels of CLC-KA protein.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of CLCNKA Antibody. The lane on the right is blocked with the CLCNKA peptide.



Western blot analysis of the lysates from HepG2 cells using CLCNKA antibody.



Western blot analysis of various cells using CLC-KA Polyclonal Antibody.