

## Anti-TNFRSF10A antibody (401-450 aa) (STJ92775)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/ELISA
Host / Source	Rabbit
Reactivity	Human/Monkey

### PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 IF 1:200-1:1000 ELISA 1:5000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Molecular Weight	Observed: 50kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION


Gene ID	<a href="#">8797</a>
Gene Symbol	<a href="#">TNFRSF10A</a>
UniProt ID	<a href="#">TR10A_HUMAN</a>
Immunogen Region	401-450 aa
Specificity	DR4 Polyclonal Antibody detects endogenous levels of DR4 protein.

### ADDITIONAL INFORMATION


Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Western blot analysis of Jurkat cells using DR4 Polyclonal Antibody



Western blot analysis of lysates from K562 cells, COS7 cells, and RAW264.7 cells, using TR10A Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of various cells using DR4 Polyclonal Antibody