

## Anti-DNM1 antibody (740-789 aa) (STJ92793)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/IF/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 IHC 1:100-1:300 ELISA 1:40000 IF 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Molecular Weight	Observed: 97kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION


Gene ID	<a href="#">1759</a>
Gene Symbol	<a href="#">DNM1</a>
UniProt ID	<a href="#">DYN1_HUMAN</a>
Immunogen Region	740-789 aa
Specificity	Dynamin I Polyclonal Antibody detects endogenous levels of Dynamin I protein.

### ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Western blot analysis of lysates from mouse brain, using Dynamin-1 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of 3T3 cells using Dynamin I Polyclonal Antibody



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using Dynamin-1 Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of various cells using Dynamin I Polyclonal Antibody