

Anti-EIF2B3 antibody (240-320 aa) (STJ92870)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 ELISA 1:5000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Molecular Weight	Observed: 50kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	8891
Gene Symbol	EIF2B3
UniProt ID	E12BG_HUMAN
Immunogen Region	240-320 aa
Specificity	eIF2B Gamma Polyclonal Antibody detects endogenous levels of eIF2B Gamma protein.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of K562 cells using eIF2B Gamma Polyclonal Antibody diluted at 1: 1000



Western blot analysis of various cells using eIF2B Gamma Polyclonal Antibody diluted at 1: 1000