

## Anti-SH3GLB2 antibody (131-180 aa) (STJ92923)

### GENERAL INFORMATION

|               |                        |
|---------------|------------------------|
| Product Type  | Primary antibodies     |
| Applications  | WB/IHC/IF/ELISA        |
| Host / Source | Rabbit                 |
| Reactivity    | Human/Mouse/Rat/Monkey |

### PRODUCT PROPERTIES


|                     |                                                                                                                       |
|---------------------|-----------------------------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                                            |
| Concentration       | 1 mg/mL                                                                                                               |
| Conjugation         | Unconjugated                                                                                                          |
| Purification        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Dilution Range      | WB 1:500-1:2000<br>IHC 1:100-1:300<br>ELISA 1:40000<br>IF 1:50-200                                                    |
| Formulation         | Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.                                               |
| Isotype             | IgG                                                                                                                   |
| Molecular Weight    | Observed: 48kD                                                                                                        |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                        |

### TARGET INFORMATION


|                  |                                                                                       |
|------------------|---------------------------------------------------------------------------------------|
| Gene ID          | <a href="#">56904</a>                                                                 |
| Gene Symbol      | <a href="#">SH3GLB2</a>                                                               |
| UniProt ID       | <a href="#">SHLB2_HUMAN</a>                                                           |
| Immunogen Region | 131-180 aa                                                                            |
| Specificity      | Endophilin B2 Polyclonal Antibody detects endogenous levels of Endophilin B2 protein. |

### ADDITIONAL INFORMATION


**Note** STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunohistochemistry analysis of paraffin-embedded human prostate carcinoma tissue, using SH3GLB2 Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from COS7 cells, using SH3GLB2 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of various cells using Endophilin B2 Polyclonal Antibody diluted at 1:1000