

## Anti-TAOK3 antibody (531-580 aa) (STJ93798)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/ELISA/IHC       |
| Host / Source | Rabbit             |
| Reactivity    | Human/Mouse/Rat    |

### PRODUCT PROPERTIES


|                     |                                                                                                                       |
|---------------------|-----------------------------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                                            |
| Concentration       | 1 mg/mL                                                                                                               |
| Conjugation         | Unconjugated                                                                                                          |
| Purification        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Dilution Range      | WB 1:500-2000<br>IHC-P 1:50-300<br>ELISA 2000-20000                                                                   |
| Formulation         | Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.                                               |
| Isotype             | IgG                                                                                                                   |
| Molecular Weight    | Observed: 110kD                                                                                                       |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                        |

### TARGET INFORMATION


|                  |                                                                   |
|------------------|-------------------------------------------------------------------|
| Gene ID          | <a href="#">51347</a>                                             |
| Gene Symbol      | <a href="#">TAOK3</a>                                             |
| UniProt ID       | <a href="#">TAOK3_HUMAN</a>                                       |
| Immunogen Region | 531-580 aa                                                        |
| Specificity      | JIK Polyclonal Antibody detects endogenous levels of JIK protein. |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of Jurkat cells using JIK Polyclonal Antibody



Western blot analysis of various cells using JIK Polyclonal Antibody