

Anti-ARHGEF12 antibody (449-498 aa) (STJ93900)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000 IHC-P 1:50-300
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Molecular Weight	Observed: 173kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	23365
Gene Symbol	ARHGEF12
UniProt ID	ARHGC_HUMAN
Immunogen Region	449-498 aa
Specificity	LARG Polyclonal Antibody detects endogenous levels of LARG protein.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Western blot analysis of RAW264.7 cells using LARG polyclonal antibody STJ93900



Western blot analysis of lysates from RAW264.7 cells, using ARHGEF12 Antibody. The lane on the right was blocked with STJ93900



Western blot analysis of the lysates from HUVEC cells using ARHGEF12 antibody. STJ93900



Immunohistochemical analysis of paraffin-embedded human liver cancer. 1, Antibody was diluted at 1:200 (4° overnight). STJ93900