

## Anti-LFNG antibody (121-170 aa) (STJ93970)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 ELISA 1:10000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Molecular Weight	Observed: 42kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION


Gene ID	<a href="#">3955</a>
Gene Symbol	<a href="#">LFNG</a>
UniProt ID	<a href="#">LFNG_HUMAN</a>
Immunogen Region	121-170 aa
Specificity	Lunatic Fringe Polyclonal Antibody detects endogenous levels of Lunatic Fringe protein.

### ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Western blot analysis of lysates from MCF-7 and HUVEC cells, using LFNG Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of 3T3 cells using Lunatic Fringe Polyclonal Antibody diluted at 1: 1000