

Anti-MERTK antibody (10-59 aa) (STJ94098)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA/IHC
Host / Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000 IHC-P 1:50-300 ELISA 2000-20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Molecular Weight	Observed: 110kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	10461
Gene Symbol	MERTK
UniProt ID	MERTK_HUMAN
Immunogen Region	10-59 aa
Specificity	MerTK Polyclonal Antibody detects endogenous levels of MerTK protein.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of HEPG2 cells using MerTK Polyclonal Antibody diluted at 1:2000



Western blot analysis of the lysates from COLO205 cells using MERTK antibody.



Western blot analysis of various cells using MerTK Polyclonal Antibody diluted at 1:2000