

## Anti-MGST2 antibody (370-450 aa) (STJ94117)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/ELISA/IHC       |
| Host / Source | Rabbit             |
| Reactivity    | Human/Rat/Mouse    |

### PRODUCT PROPERTIES


|                     |                                                                                                                       |
|---------------------|-----------------------------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                                            |
| Concentration       | 1 mg/mL                                                                                                               |
| Conjugation         | Unconjugated                                                                                                          |
| Purification        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Dilution Range      | WB 1:500-2000<br>IHC-P 1:50-300<br>ELISA 2000-20000                                                                   |
| Formulation         | Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.                                               |
| Isotype             | IgG                                                                                                                   |
| Molecular Weight    | Observed: 17kD                                                                                                        |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                        |

### TARGET INFORMATION


|                  |                                                                       |
|------------------|-----------------------------------------------------------------------|
| Gene ID          | <a href="#">4258</a>                                                  |
| Gene Symbol      | <a href="#">MGST2</a>                                                 |
| UniProt ID       | <a href="#">MGST2_HUMAN</a>                                           |
| Immunogen Region | 370-450 aa                                                            |
| Specificity      | MGST2 Polyclonal Antibody detects endogenous levels of MGST2 protein. |

### ADDITIONAL INFORMATION

**Note** STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of 293T cells using MGST2 polyclonal antibody diluted at 1:2000 STJ94117



Immunohistochemical analysis of paraffin-embedded human liver cancer. 1. Antibody was diluted at 1:200 (4°C overnight). STJ94117