

Anti-SLC22A3 antibody (275-324) (STJ94597)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 ELISA 1:10000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Molecular Weight	Observed: 61kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	6581
Gene Symbol	SLC22A3
UniProt ID	S22A3_HUMAN
Immunogen Region	275-324
Specificity	OCT3 Polyclonal Antibody detects endogenous levels of 40819 protein.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



CPSF3 Promotes Pre-mRNA Splicing and Prevents circRNA Cyclization in Hepatocellular Carcinoma. Cancers, Qubo Zhu, WB, STJ94597



Western Blot analysis of various cells using OCT3 Polyclonal Antibody STJ94597



Western blot analysis of lysates from HUVEC, COLO, and HT-29 cells, using SLC22A3 Antibody. The lane on the right is STJ94597