

Anti-PARD3 antibody (1141-1190 aa) (STJ94955)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/IF/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 IHC 1:100-1:300 IF 1:200-1:1000 ELISA 1:20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Molecular Weight	Observed: 151kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	56288
Gene Symbol	PARD3
UniProt ID	PARD3_HUMAN
Immunogen Region	1141-1190 aa
Specificity	PARD3A Polyclonal Antibody detects endogenous levels of PARD3A protein.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Western Blot analysis of various cell lines using PARD3A Polyclonal Antibody STJ94955



Western blot analysis of lysates from COLO205 cells, using PARD3 Antibody. The lane on the right is blocked with the STJ94955



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200 (4° overnight). 2, STJ94955



Immunofluorescence analysis of HepG2 cells, using PARD3 Antibody. The picture on the right is blocked with the STJ94955