

Anti-PIK3R2 antibody (409-458 aa) (STJ95080)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 ELISA 1:10000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Molecular Weight	Observed: 85kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	5296
Gene Symbol	PIK3R2
UniProt ID	P85B_HUMAN
Immunogen Region	409-458 aa
Specificity	PI 3-kinase p85 Beta Polyclonal Antibody detects endogenous levels of PI 3-kinase p85 Beta protein.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of 453 cells using PI 3-kinase p85 Beta Polyclonal Antibody diluted at 1: 500



Western blot analysis of lysate from Jurkat cells, using PI 3-kinase p85 Beta antibody.



Western blot analysis of various cells using PI 3-kinase p85 Beta Polyclonal Antibody diluted at 1: 500