

Anti-PPP2R5D antibody (544-593 aa) (STJ95194)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/IF/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 IHC 1:100-1:300 ELISA 1:10000 IF 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Molecular Weight	Observed: 70kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	5528
Gene Symbol	PPP2R5D
UniProt ID	2A5D_HUMAN
Immunogen Region	544-593 aa
Specificity	PP2A-B56-Delta Polyclonal Antibody detects endogenous levels of PP2A-B56-Delta protein.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Western Blot analysis of 293T cells using PP2A-B56-Δ Polyclonal Antibody diluted at 1:1000 STJ95194



Western blot analysis of lysates from Jurkat and HUVEC cells, using PPP2R5D Antibody. The lane on the right was blocked STJ95194



Immunohistochemical analysis of paraffin-embedded human brain tumor. 1, Antibody was diluted at 1:200 (4° overnight). STJ95194



Immunofluorescence analysis of SiHa cells. 1, Primary antibody was diluted at 1:100 (4°C overnight). 2, Goat Anti Rabbit STJ95194