

## Anti-TMBIM1 antibody (251-300 aa) (STJ95403)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/IHC/IF/ELISA    |
| Host / Source | Rabbit             |
| Reactivity    | Human/Mouse/Rat    |

### PRODUCT PROPERTIES


|                     |                                                                                                                       |
|---------------------|-----------------------------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                                            |
| Concentration       | 1 mg/mL                                                                                                               |
| Conjugation         | Unconjugated                                                                                                          |
| Purification        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Dilution Range      | WB 1:500-1:2000<br>IHC 1:100-1:300<br>ELISA 1:40000<br>IF 1:50-200                                                    |
| Formulation         | Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.                                               |
| Isotype             | IgG                                                                                                                   |
| Molecular Weight    | Observed: 35kD                                                                                                        |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                        |

### TARGET INFORMATION


|                  |                                                                       |
|------------------|-----------------------------------------------------------------------|
| Gene ID          | <a href="#">64114</a>                                                 |
| Gene Symbol      | <a href="#">TMBIM1</a>                                                |
| UniProt ID       | <a href="#">LFG3_HUMAN</a>                                            |
| Immunogen Region | 251-300 aa                                                            |
| Specificity      | RECS1 Polyclonal Antibody detects endogenous levels of RECS1 protein. |

### ADDITIONAL INFORMATION


**Note** STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Western Blot analysis of various cells using RECS1 Polyclonal Antibody diluted at 1:1000 STJ95403



Western blot analysis of lysates from HUVEC cells, using TMBIM1 Antibody. The lane on the right was blocked with the STJ95403



Western blot analysis of the lysates from K562 cells using TMBIM1 antibody. STJ95403



Immunohistochemical analysis of paraffin-embedded human colon cancer. 1, Tris-EDTA, pH 9.0 was used for antigen STJ95403