

Anti-RPTN antibody (648-697 aa) (STJ95417)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC
Host / Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000 IHC-P 1:50-300
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Molecular Weight	Observed: 90kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	126638
Gene Symbol	RPTN
UniProt ID	RPTN_HUMAN
Immunogen Region	648-697 aa
Specificity	Repetin Polyclonal Antibody detects endogenous levels of Repetin protein.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Western Blot analysis of HeLa cells using Repetin Polyclonal Antibody STJ95417



Western blot analysis of RPTN Antibody. The lane on the right is blocked with the RPTN peptide. STJ95417



Western blot analysis of the lysates from HepG2 cells using RPTN antibody. STJ95417



Immunohistochemical analysis of paraffin-embedded human breast cancer. 1. Antibody was diluted at 1:200 (4° overnight). STJ95417