

Anti-SDHA antibody (551-600 aa) (STJ95591)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000 IHC-P 1:50-300
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Molecular Weight	Observed: 70kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	6389
Gene Symbol	SDHA
UniProt ID	SDHA_HUMAN
Immunogen Region	551-600 aa
Specificity	SDHA Polyclonal Antibody detects endogenous levels of SDHA protein.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of HepG2 cells using SDHA Polyclonal Antibody diluted at 1:2000



Western blot analysis of lysate from HeLa cells, using SDHA antibody.



Western blot analysis of various cells using SDHA Polyclonal Antibody diluted at 1:2000