

## Anti-HTR1E antibody (101-150 aa) (STJ95759)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IF/ELISA/IHC
Host / Source	Rabbit
Reactivity	Human

### PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	IF ICC 1:50-300 ELISA 1:2000-20000 IHC-P 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Molecular Weight	Observed: 42kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION


Gene ID	<a href="#">3354</a>
Gene Symbol	<a href="#">HTR1E</a>
UniProt ID	<a href="#">5HT1E_HUMAN</a>
Immunogen Region	101-150 aa
Specificity	SR-1E Polyclonal Antibody detects endogenous levels of SR-1E protein.

### ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western Blot analysis of various cells using SR-1E Polyclonal Antibody STJ95759



Western blot analysis of lysates from HeLa cells, using 5-HT-1E Antibody. The lane on the right was blocked with the STJ95759



Immunohistochemical analysis of paraffin-embedded human squamous cell carcinoma of the lung. 1, Antibody was diluted at STJ95759