

Anti-TBX3 antibody (301-350 aa) (STJ95936)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 ELISA 1:5000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Molecular Weight	Observed: 79kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	6926
Gene Symbol	TBX3
UniProt ID	TBX3_HUMAN
Immunogen Region	301-350 aa
Specificity	TBX3 Polyclonal Antibody detects endogenous levels of TBX3 protein.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of mouse-kidney, mouse-brain, KB, 293T lysis using TBX3 antibody. Polyclonal Antibody diluted at 1:1000



Western blot analysis of KB cells using TBX3 Polyclonal Antibody diluted at 1:1000