

## Anti-ZNF575 antibody (10-59 aa) (STJ96337)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/ELISA           |
| Host / Source | Rabbit             |
| Reactivity    | Human/Rat/Mouse    |

### PRODUCT PROPERTIES


|                     |                                                                                                                       |
|---------------------|-----------------------------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                                            |
| Concentration       | 1 mg/mL                                                                                                               |
| Conjugation         | Unconjugated                                                                                                          |
| Purification        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Dilution Range      | WB 1:500-1:2000<br>ELISA 1:5000                                                                                       |
| Formulation         | Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.                                               |
| Isotype             | IgG                                                                                                                   |
| Molecular Weight    | Observed: 26kD                                                                                                        |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                        |

### TARGET INFORMATION


|                  |                                                                         |
|------------------|-------------------------------------------------------------------------|
| Gene ID          | <a href="#">284346</a>                                                  |
| Gene Symbol      | <a href="#">ZNF575</a>                                                  |
| UniProt ID       | <a href="#">ZNF575_HUMAN</a>                                            |
| Immunogen Region | 10-59 aa                                                                |
| Specificity      | ZNF575 Polyclonal Antibody detects endogenous levels of ZNF575 protein. |

### ADDITIONAL INFORMATION

**Note** STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of the lysates from HeLa cells using ZNF575 antibody.



Western blot analysis of ZNF575 Antibody. The lane on the right is blocked with the ZNF575 peptide.