

Anti-ILF3 antibody (302-351 aa) (STJ96401)

GENERAL INFORMATION

	Product Type	Primary antibodies
	Applications	WB/ELISA
	Host / Source	Rabbit
	Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 ELISA 1:20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Molecular Weight	Observed: 95kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	3609
Gene Symbol	ILF3
UniProt ID	ILF3_HUMAN
Immunogen Region	302-351 aa
Specificity	NF90 Polyclonal Antibody detects endogenous levels of NF90 protein.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Western Blot (WB): Zhuang, Junlong, et al.
"TGFβ1 promotes gemcitabine resistance through
regulating the STJ96401



Western Blot analysis of extracts from A549, K562, NIH-3T3 cells, using NF90 Polyclonal Antibody. Secondary antibody STJ96401



Western blot analysis of lysates from HeLa cells, using NF90 antibody. STJ96401



Immunofluorescence analysis of MCF7 cells. 1, Primary antibody was diluted at 1:100 (4°C overnight). 2, Goat Anti Rabbit STJ96401