

## Anti-IL5RA antibody (41-90 aa) (STJ96579)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/IF/ELISA
Host / Source	Rabbit
Reactivity	Human/Rat/Mouse

### PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 IHC-P 1:100-300 ELISA 1:20000 IF 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Molecular Weight	Observed: 47kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION


Gene ID	<a href="#">3568</a>
Gene Symbol	<a href="#">IL5RA</a>
UniProt ID	<a href="#">IL5RA_HUMAN</a>
Immunogen Region	41-90 aa
Specificity	IL-5R Alpha Polyclonal Antibody detects endogenous levels of IL-5R Alpha protein.

### ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Western blot analysis of lysate from K562 cells, using IL5RA Antibody.



Western blot analysis of K562 cells using IL-5R Alpha Polyclonal Antibody.. Secondary antibody was diluted at 1:20000