

## Anti-CGRRF1 antibody (283-332 aa) (STJ96731)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/ELISA           |
| Host / Source | Rabbit             |
| Reactivity    | Human/Rat/Mouse    |

### PRODUCT PROPERTIES


|                     |                                                                                                                       |
|---------------------|-----------------------------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                                            |
| Concentration       | 1 mg/mL                                                                                                               |
| Conjugation         | Unconjugated                                                                                                          |
| Purification        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Dilution Range      | WB 1:500-1:2000<br>ELISA 1:20000                                                                                      |
| Formulation         | Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.                                               |
| Isotype             | IgG                                                                                                                   |
| Molecular Weight    | Observed: 38kD                                                                                                        |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                        |

### TARGET INFORMATION


|                  |                                                                       |
|------------------|-----------------------------------------------------------------------|
| Gene ID          | <a href="#">10668</a>                                                 |
| Gene Symbol      | <a href="#">CGRRF1</a>                                                |
| UniProt ID       | <a href="#">CGRF1_HUMAN</a>                                           |
| Immunogen Region | 283-332 aa                                                            |
| Specificity      | CGR19 Polyclonal Antibody detects endogenous levels of CGR19 protein. |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of lysate from A549 cells, using CGRRF1 Antibody.



Western blot analysis of A549 cells using CGR19 Polyclonal Antibody.. Secondary antibody was diluted at 1:20000