

## Anti-LAIR1 antibody (21-70 aa) (STJ97329)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Rat/Mouse

### PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 ELISA 1:10000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Molecular Weight	Observed: 32kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION


Gene ID	<a href="#">3903</a>
Gene Symbol	<a href="#">LAIR1</a>
UniProt ID	<a href="#">LAIR1_HUMAN</a>
Immunogen Region	21-70 aa
Specificity	CD305 Polyclonal Antibody detects endogenous levels of CD305 protein.

### ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Western blot analysis of KB cells using CD305 Polyclonal Antibody. Secondary antibody was diluted at 1:20000



Western blot analysis of KB, NIH-3T3 cells using CD305 Polyclonal Antibody. Secondary antibody was diluted at 1:20000