

Anti-PTGFRN antibody (601-650 aa) (STJ97332)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/IF/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 IHC-P 1:100-300 ELISA 1:20000 IF 1:100-300
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Molecular Weight	Observed: 98kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	5738
Gene Symbol	PTGFRN
UniProt ID	FPRP_HUMAN
Immunogen Region	601-650 aa
Specificity	CD315 Polyclonal Antibody detects endogenous levels of CD315 protein.

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western Blot analysis of 293 cells using CD315 Polyclonal Antibody. Secondary antibody was diluted STJ97332



Immunohistochemical analysis of paraffin-embedded human liver, antibody was diluted at 1:200 STJ97332



Immunofluorescence analysis of A549. 1, primary antibody was diluted at 1:200 (4°C overnight). 2, Goat Anti Rabbit IgG STJ97332