

Anti-FABP4 antibody (61-121 aa) [5H11A4E11] (STJ97815)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 ELISA 1:10000
Formulation	Liquid in PBS containing 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG1
Molecular Weight	Calculated: 15kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	2167
Gene Symbol	FABP4
UniProt ID	FABP4_HUMAN
Immunogen Region	61-121 aa
Specificity	A-FABP Monoclonal Antibody detects endogenous levels of A-FABP protein.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis using A-FABP monoclonal antibody against truncated Trx-FABP4 recombinant protein (1). STJ97815



Western Blot (WB): A Novel Role in Duck Myoblasts: The Transcription Factor Retinoid X Receptor Alpha (RXRA) Inhibits STJ97815