

## Anti-APOA5 antibody [1G5G9 (c) /2G1H11G9B3] (STJ97836)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES


Clonality	Monoclonal
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 ELISA 1:10000
Formulation	Liquid in PBS containing 0.03% Sodium Azide.
Isotype	IgG1
Molecular Weight	Calculated: 41kD
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">116519</a>
Gene Symbol	<a href="#">APOA5</a>
UniProt ID	<a href="#">APOA5_HUMAN</a>
Specificity	ApoA-V (ab) Monoclonal Antibody detects endogenous levels of ApoA-V (ab) protein.

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis using ApoA-V (ab) monoclonal antibody against human serum (1) and ApoA5 recombinant protein (2).